

where:  $m$ =moisture content of the sample.

**§ 436.301 Thin layer chromatography identity test for carbenicillin indanyl.**

Using the sample solution prepared as described in the section for the antibiotic drug to be tested, proceed as described in paragraphs (a), (b), (c), and (d) of this section.

(a) *Equipment*—(1) *Chromatography tank*. A rectangular tank, approximately  $9 \times 9 \times 3.5$  inches lined with Whatman's 3MM chromatographic paper (0.3 millimeters) or equivalent.

(2) *Iodine vapor chamber*. A rectangular tank approximately  $9 \times 9 \times 3.5$  inches, with a suitable cover, containing iodine crystals.

(3) *Plates*. Use  $20 \times 20$  centimeters thin layer chromatography plates coated with silica gel G or equivalent to a thickness of 250 microns.

(b) *Reagents*—(1) *Extraction solvent*. Mix ethyl acetate, acetone, pyridine, water, and acetic acid in volumetric proportions of 100:200:25:75:1.5 respectively.

(2) *Developing solvent*. Mix ethyl acetate, acetone, pyridine, water, and acetic acid in volumetric proportions of 300:400:25:75:2 respectively.

(3) *Ferric chloride-potassium ferricyanide reagent*. Immediately before use, mix 100 milliliters of a 1 percent ferric chloride solution in 1 percent hydrochloric acid with 100 milliliters of a 1 percent potassium ferricyanide solution and 75 milliliters of methanol.

(c) *Preparation of working standard solution*. Weigh an amount of the carbenicillin indanyl working standard equivalent to approximately 10 milligrams of carbenicillin into a 50-milliliter Erlenmeyer flask. Dissolve the material in sufficient extraction solvent to make a solution containing 1 milligram carbenicillin per milliliter.

(d) *Procedure*. Pour developing solvent into the bottom of the chromatography tank. Cover and seal the tank. Allow it to equilibrate for 1 hour. Prepare a plate as follows: On a line 2 centimeters from the base of the silica gel plate, and at intervals of 2 centimeters, spot 10 microliters of the standard solution and the sample solution. The plate should be air dried for 30 minutes. Place the plate into the

chromatography tank. Allow the solvent front to travel about 15 centimeters from the starting line and then remove the plate from the tank. Heat the plate for 30 minutes at  $80^\circ\text{C}$ . in a circulating air oven and then allow the plate to cool to room temperature. Place the plate in the iodine vapor chamber for about 30 seconds, remove the plate and spray it with the ferric chloride-potassium ferricyanide reagent. Carbenicillin indanyl appears as a blue spot on a yellow-green background at an  $R_f$  of about 0.5. The test is satisfactory if the sample compares qualitatively with the standard.

[39 FR 18944, May 30, 1974, as amended at 41 FR 18509, May 5, 1976]

**§ 436.302 Clindamycin vapor phase chromatography.**

(a) *Equipment*. Gas chromatograph equipped with a flame ionization detector: Barber-Colman 5,000 or equivalent.

(b) *Reagents*. (1) Pyridine, reagent grade, dried over sodium sulfate.

(2) Chloroform, reagent grade.

(3) Acetic anhydride, reagent grade, used as acetylating agent.

(4) Internal standard: Prepare a solution containing 3 milligrams of cholestane per milliliter in pyridine.

(c) *Typical conditions*. (1) Column: 4 feet  $\times$  4 millimeters ID, glass, with 1 percent SE-30 on Diatoport S (60/80 mesh), or equivalent.

(2) Temperatures: Column  $200^\circ\text{C}$ .; detector  $215^\circ\text{C}$ .; injection port, ambient temperature.

(3) Carrier gas: Helium approximately 120 milliliters per minute.

(4) Detector: Hydrogen flame—hydrogen at 120 pounds per square inch, air at 40 pounds per square inch.

(5) Sensitivity: 1,000; attenuation, 2 for clindamycin, 1 for internal standard:  $2 \times 10^{-8}$  amperes.

(d) *Preparation of clindamycin sample and working standard solutions*. Accurately weigh approximately 15 milligrams of sample or working standard into a glass-stoppered conical 15-milliliter centrifuge tube. Add 1.0 milliliter of chloroform, 1.0 milliliter of internal standard solution, and 0.6 milliliter of acetic anhydride. Agitate the tubes to insure dissolution of the sample and

complete mixing of the liquids. Proceed as directed in paragraph (e) of this section.

(e) *Procedure.* Cover the top of each centrifuge tube with a plastic cap. Punch a small hole in the top of each cap to allow vapor to escape. Place the tubes in a 100° C. drying oven for 2.5 hours. Remove the tubes from the oven and allow to cool. Take the plastic cap from each tube and replace with the glass stopper. Centrifuge 10–15 minutes at 2,000–2,500 r.p.m. to separate the white solid from the liquid in the tube. Inject 0.5 microliter of the clear liquid into the gas chromatograph. Use the conditions and materials listed in paragraphs (a), (b), and (c) of this section. The conditions should be adequate to maintain a stable baseline and provide at least 60 percent deflection of the recorder scale by the clindamycin peak. The resolution of the peaks should be complete. The elution order is: Internal standard, clindamycin, and epiclindamycin (if present). Calculate the clindamycin content as directed in paragraph (f) of this section.

(f) *Calculations.* Calculate the clindamycin content of the sample as follows:

$$\text{Micrograms of clindamycin per milligram} = \frac{R_u \times W_s \times f}{R_s \times W_u}$$

where:

$R_u$  = Area of the clindamycin sample peak (at a retention time equal to that observed for the clindamycin standard)/Area of internal standard peak;

$R_s$  = Area of the clindamycin standard peak/Area of internal standard peak;

$W_s$  = Weight of the clindamycin working standard in milligrams;

$W_u$  = Weight of the sample in milligrams;

$f$  = Potency of the clindamycin working standard in micrograms per milligram.

**§ 436.303 Clindamycin content of clindamycin palmitate hydrochloride by vapor phase chromatography.**

(a) *Equipment.* Gas chromatograph equipped with a flame ionization detector: Hewlett-Packard 7606<sup>4</sup> or equivalent.

(b) *Reagents.* (1) Acetic anhydride, reagent grade.

(2) Pyridine, reagent grade.

(3) Chloroform, reagent grade.

(4) Internal standard: Prepare a solution containing 5 milligrams of cholesteryl benzoate per milliliter in chloroform.

(c) *Typical conditions.* (1) Column: 6 feet × 2 millimeters ID, glass, with 1 percent UC-W98 on Chromosorb WHP (80/100 mesh) or equivalent.

(2) Temperatures: Column 275° C.; detector 290° C.; injection port 280° C.

(3) Carrier gas: Helium approximately 60 milliliters per minute.

(4) Detector: Hydrogen flame ionization—hydrogen at 12 pounds per square inch, air at 32 pounds per square inch.

(5) Sensitivity: 1,000; attenuation, 16;  $1 \times 10^{-9}$  amperes.

(d) *Preparation of clindamycin palmitate hydrochloride sample and working standard solutions.* Accurately weigh approximately 15 milligrams of both the sample and the working standard into separate glass-stoppered, conical 15-milliliter centrifuge tubes. Add 1.0 milliliter of internal standard solution, 1.0 milliliter of pyridine, and 0.5 milliliter of acetic anhydride to each tube. Agitate the tubes to insure dissolution and complete mixing of the liquids. Proceed as directed in paragraph (e) of this section.

(e) *Procedure.* Cover the top of each centrifuge tube with a plastic cap. Punch a small hole in the top of each cap to allow vapor to escape. Place the tubes in a 100° C. drying oven for 2.5 hours. Remove the tubes from the oven and allow to cool. Take the plastic cap from each tube and replace with the glass stopper. Centrifuge 10–15 minutes at 2,000–2,500 r.p.m. to separate the white solid from the liquid in the tube. Inject 1 microliter of the clear liquid into the gas chromatograph. Use the conditions and materials listed in paragraphs (a), (b), and (c) of this section. The conditions should be adequate to maintain a stable baseline and provide at least 40 percent deflection of the recorder scale by the clindamycin palmitate peak. The resolution of the peaks should be complete. The internal standard will be eluted before the clindamycin palmitate. Calculate the clindamycin content as directed in paragraph (f) of this section.

<sup>4</sup> Available from: Hewlett Packard Co., P.O. Box 301, Loveland, CO 80537.